**ČIUŽINYS PRAGULŲ PROFILAKTIKAI
MATRACIS IZGULĖJUMU PROFILAKSEI
KOMPRESSORIGA MADRATS LAMATISTE PROFŪLAKTIKAKS
BUBBLE ANTIDECUBITUS MATTRESS
МАТРАС С КОМПРЕССОРОМ ДЛЯ ПРОФИЛАКТИКИ ПРОЛЕЖНЕЙ**

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PREFACE

Thank you for the confidence shown in Kid-Man. Every product is factory checked and left us in faultless condition, We wish you much success and we hope Kid-Man product will make a positive impact to your daily living.

INDICATIONS

A set, pump and mattress, is developed for bed sore proflaxis and treatment. Device is intended to be used both in domestic and hospital environment to prevent bedridden people of developing bed sores. Using this device in improper way may cause malfunction of device or damage to product or user.

SYMBOLOLOGY

Caution

Do not bleach

Medical Device Directive 93/42/EEC

Electrical Protection Type B

Tumble dry with low heat

Class II Equipment (Double Insulated)

Date of manufacture

Do not iron

Refer to manual before use

Compatible with any dry cleaning methods

The marking of electrical and electronics devices according to Directive 2002/96/EC. The device accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances of Regulations for disposal.

No more than 80° C wash

1. Do not use this product or any available optional equipment without first completely reading and understanding this instruction manual. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, dealer or technical personnel before attempting to use this equipment, otherwise injury or damage may occur.

2. Advise with physician or therapist to determine the correct adjustment and the correct use of the device.

3. Keep the packed kit away from heat sources.

4. The device use limit is defined by the wear of the parts.

5. Do not allow children to play on or operate the pump.

6. Do not use device in other way than it is intended to be used

7. UAB Kasko Group is not responsible for damage occurred due to improper use of device or unauthorized repair.

8. Manufacturer has right to change the information in this user manual before prior notice

- Name of products
- Model
- Maximum load
- Manufacturer
- Contact details
- Serial number

Pic. No. 25

Pic. No. 26


PARTS OF THE PRODUCT

02-M8: mattress and pump.

NOTE: Pictures are only for the demonstration and do not refer to a specific product.

CONTROL AND FEATURES

Alternating pressure system

Pic. No. 27

MAIN CONTROL PANEL

1. Power switch (1)

Press this button to turn on the pump.

2. Alternating/Static (2)

Alternating: normal alternating function.

Static: To keep mattress always inflating


3. Adjustable valve (3)

Adjust the valve to increase or decrease for a softer or firmer setting between 1-8 range.


Control panel

Pic. No. 28

Front:


1. Control panel
2. Female quick connector
3. Male quick connector

Pic. No. 29

Back:

4. Filter
5. Main filter
6. Protective cap of filter
7. Hook

Pic. No. 30


1. Cover
2. Fixing belt
3. Cells (can be taken off)
4. Male quick connector
5. Air tubing

SET UP

1. Place the mattress on a common mattress, with the air hose connectors oriented to the footboard.

NOTE: lean the mattress on a 5cm mattress.

2. Hang the pump onto the bedrail (footboard side), or place it on a flat surface.

Pic. No. 31

3. Connect air hose connectors from air mattress to the pump unit. Ensure the air hoses are not kinked or tucked under mattress.

4. Plug the power cord into the electrical outlet.

5. Turn the switch ON.

6. Check the correct assembly, paying particular attention to air hoses connector before starting using.

Pic. No. 32

Make sure the pump unit is suitable for the local power voltage.
Do not position the equipment so it difficult to disconnecting the plug.

USING WARNINGS

- Do not keep pump in place where it would be difficult to turn it off easily!
- Please unplug the control unit, if do not need to use for a long time!
- Avoid the cable, plug and main control unit contact with water or other fluid!
- Avoid cables damage when moved the sickbed!

- Ensure that cables are safe and tidy and not caught in the bed mechanism or trailing on the floor
- Avoid sharp objects or strong acid and alkali chemical agents from mattress.
- Avoid use the system with connecting power without any people watch
- Ensure the mattress, cover, cables and power outlet out of fire, heat or strong acid and alkali chemical agents
- Please plug the control unit into the ground power supply
- Do not modify the pump and mattress without authorization of the manufacturer and do not stuff other objects into the crack of control unit or mattress. Please don't use the parts which not belong to Kasko Group, UAB.

DO NOT USE DEVICE IF:

- The system sounds abnormal or vibrate intensively
- Control unit drop on the floor and damaged.
- Control unit drop into the water or other fluid.

USE OF THE PRODUCT

Patient weight	Position on the pressure knob
<50 kg	2-3
50 - 70 kg	3-4
70 - 90 kg	4-5
90 - 110 kg	5-6
110 - 130 kg	6-7
130 - 150 kg	8

NOTE: In case the pressure is constantly low, check for any holes (on the interchangeable elements or on the connection tubes). If it is necessary substitute the damaged connection tubes or the interchangeable elements, otherwise contact your dealer. If a rapid deflation is necessary unplug the connection tube from the mattress.

STORAGE

If mattress is not being used for a long period, please prepare it for storage.

- 1) Switch off power from main and unplug, separate the pump with Mattress and put the mattress flat.
- 2) Fold the mattress in two from mattress foot to head.
- 3) Allow time for air to escape from mattress, then reduce the space to storage.
- 4) Closed the cover the protective cap on the quick connector.

MAINTENANCE

The Pump & Mattress devices by UAB Kasko Group are checked carefully and supplied with CE mark, once launched on the market.

For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every year. In case of repair, use only original spare parts and accessories.

CLEANING AND DISINFECTION

WARNING!

Turn pump off before starting cleaning procedure. Never use acids, alkalis or solvents such as acetone.

1. Clean the mattress with mild soap.
2. Do not immerse or soak pump unit.
3. After cleaning, dry the mattress without direct sunlight exposure.
4. Check the connection tubes are not obstructed by dust or else. In this case use a soft swab to clean the tubes.

DISPOSAL

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the kit in the appropriate disposal areas for recycling.

BATTERIES (DIRECTIVE 2002/96/EC)

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the product is marked with crossed mobile waste bin.

TROUBLESHOOTING

If your questions can't be answered with below information, please contact your local agent directly. They might require a technician to take care of the problem.

NO LIGHTS ON UNIT	<ul style="list-style-type: none"> • Check control unit is connect to the mains power supply and that the unit is switched on. • Check quick blow fuse on the control unit or mains plug fuse.
LOW PRESSURE ALERT	<ul style="list-style-type: none"> • Check main power to the control unit and that the plug is switched. • Check mattress air supply connector is connected to the pump unit correctly. Check CPR connector is connected properly. • Under the mattress cover check all "T" & "L" connectors are connected to the air tubes. • Check the connector tubes for kinks, crimps or damage
PATIENT'S BODY IS SAGGING IN THE MIDDLE/ BOTTOMING OUT	<p>Increase the value of the comfort the setting. Wait a couple of minutes so pressure can stabilize before making another change.</p>
CONTROL LOCK UP OR "FREEZE"	<p>Turn the control unit off and unplug the pump. Wait a few minutes and plug the control unit back in. Turn on the control unit.</p>

ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturers declaration action-electromagnetic emission

The device is intended for use in the electromagnetic environment specified below.
The customer or the user of the device should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions IEC 61000-3-3	Compliance	

Guidance and manufacturers declaration - electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IFC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input / output lines	+/- 2 kV for power supply lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1 kV line(s) to line(s) +/- 2 kV line(s) to earth	+/- 1 kV differential mode +/- 2 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	<5% UT (>95% dip in UT) 0,5 cycle 40% UT (60% dip in UT) 5 cycles 70% UT (30% dip in UT) 25 cycles <5% UT (>95% dip in UT) 5 s	<5% UT (>95% dip in UT) 0,5 cycle 40% UT (60% dip in UT) 5 cycles 70% UT (30% dip in UT) 25 cycles <5% UT (>95% dip in UT) 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptable power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The power frequency magnetic fields should be at levels characteristic of a typical commercial or hospital environment.

NOTE UT is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturers declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below.

The customer or the user of the device should assure that is used in such an environment.

Immunity test	IFC 60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz – 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance: $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$ 80 MHz to 800 MHz $d=2,3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^A, should be less than the compliance level in each frequency range^B. Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2,5 GHz	3 V/m	

NOTE! At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE! These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast can not be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeded the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorientating or relocating the device.

B Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz – 80 MHz $d=1,2 \times P$	80 MHz – 800 MHz $d=1,2 \times P$	800 MHz – 2,5 GHz $d=1,2 \times P$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE! At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE! These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

SPECIFICATIONS

KOMPRESORIUS	
Power supply	120V/60Hz, 230V/50Hz
Power	6,0 W
Fuse rating	500 mA / 250 V
Air pressure	30 mmhg ~ 60 mmhg
Dimensions	23,5 cm x 26 cm x 11,9 cm
Cycle time	10 min
Weight	2,3 Kg
Environment temperature	Operation +10°C / +40°C Storage -15°C/+50°C Shipping -15°C/+70°C
Environment humidity	Operation 10%-90% Storage 10%-90% Shipping 10%-90%
Athmospheric pressure	Operation 700-1013,25 hPa

MATTRESS	
Dimensions (without flaps)	86 cm x 196 cm x 12,7 cm
Cell height (inflated)	12,7 cm
Air cell material	PVC
Cover material	Nylon
Elements number	16
Maximum patient weight	150 kg
Weight	5,5 kg

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• We issue a warranty of 24 months from date of purchase. Any unauthorized modifications carried out without approval of the manufacturer will make the liability null and void. This includes exceeding the maximum user loads and using product in improper way. Warranty is only valid with fully filled warranty card. UAB „Kasko Group“ after getting claimed item from purchaser within period of 4 weeks decides whether issue is treated under warranty. If issue is not covered under warranty and customer agrees, UAB „Kasko Group“, if possible, fixes issue. All costs, including costs of transportation, should be covered by customer. Warranty repair does not extend warranty.

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